

hrmcom.txt

```
////////////////////////////////////
//
////////
//
//  SYSTEM:          Polar HRMCom
//
//  UNIT FILE NAME:  HRMCOM.H
//
//  MODULE:          HRMCOM.DLL
//
//  AUTHORS:         MEr / Polar Electro Oy
//
//  VERSION:         1.2
//
//  DATE:            20.03.2001
//
//  ABSTRACT:        Main header file for Polar HRMCOM.DLL function
library file.
//
//  REMARKS:
//
//  COPYRIGHT (C) 2001 BY POLAR ELECTRO OY
//
////////////////////////////////////
//
////////
```

□

```
////////////////////////////////////
//
////////
//
//  BOOLEAN VARIABLES
//  -----
//
//  Function library uses boolean variables as following:
//
//          TRUE      equals  1
//          FALSE     equals  0
//
////////////////////////////////////
//
////////

////////////////////////////////////
//
////////
```


hrmcom.txt

```

x          x          x          x4          x4
//      Exercise Sets          x          x          x

x          x          x1         x4          x4
//      User Settings          x          x          x

x2         x2         x2         x4          x4
//      Reminders

x          x          x
//      Exercise Profiles

x
//      Monitor Bitmaps

x          x          x
//      Bikes

x          x3
//
//      x          = feature available
//      x1         = Only one exercise set available
//      x2         = User settings extended with user name string
//      x3         = Also power output settings
//      x4         = Education models' features are limited, see function
definitions
//
//      For more details about feature difference, refer to each
function call definitions
//      and HR monitor user's manuals. Another good hint is also to us
e
Polar Precision
//      Performance SW 3.0 and it's HR Monitor Connection. This softwa
re
products utilizes
//      HRMCom.dll function library.
//
////////////////////////////////////
//
////////
//
////////////////////////////////////
//
////////
//
// POLAR CHARACTER SET

```

hrmcom.txt

```
// -----
//
//      The following characters are valid at text strings in settings
//      :
//      - Capital letters:      ABCDEFGHIJKLMNOPQRSTUVWXYZ
//      - Small letters:       abcdefghijklmnopqrstuvwxyz
//      - Numbers:              0123456789
//      - Special chars:       -%/( ) * + . : ? and space
//
//      Unrecognized characters will be converted automatically to
spaces.
//      The text strings have to be ended by zero character (NULL).
//      Strings can be checked before sending by using function
fnHRMCom_CheckPolarCharString
//
__declspec (dllexport) BOOL CALLBACK fnHRMCom_CheckPolarCharString
(LPTSTR);
//
////////////////////////////////////
//
////////
//
////////////////////////////////////
//
////////
//
//      DATE FORMAT
//      -----
//
//      Date values are processed in yyyyymmdd format:
//
//          yyyy      year      4 digits
//          mm         month    2 digits
//          dd         day       2 digits
//
//      For example:   August 2nd 2000      => 20000802
//                      December 24th 2003  =>
20031224
//
//      NOTE: Leading zero with days and months is always obligatory.
//
////////////////////////////////////
//
////////
//
////////////////////////////////////
//
////////
```

hrmcom.txt

```
//
//      ERROR CHECKING
//      -----
//
//      All input parameters will be checked before sending to heart
//      rate monitor.
//      If any erratic values are determined, function call returns
//      FALSE and does not
//      continue sending data to monitor. The latest error code can be
//
//      checked by function:
//      iError = fnHRMCom_GetErrorCode ();
//
//      __declspec (dllexport) int CALLBACK fnHRMCom_GetErrorCode
//
//      (void);
//
//      ## UNDER CONSTRUCTION ## //
//
//      //////////////////////////////////////
//
//      //////////////////////////////////////
//
//      //////////////////////////////////////
//
//      //////////////////////////////////////
//
//      POLAR HR MONITOR TYPES
//      -----
//
//      #define      HRM_S210
//      8              // Polar S210
//      #define      HRM_S410
//      9              // Polar S410
//      #define      HRM_S510
//      10             // Polar S510
//      #define      HRM_S610
//      11             // Polar S610
//      #define      HRM_S710
//      12             // Polar S710
//      #define      HRM_S810
//      13             // Polar S810
//      #define      HRM_E200
//      14             // Polar E200 Education HRM
//      #define      HRM_E600
//      15             // Polar E600 Education HRM
//
//      NOTE: Education HR monitors are regionally available.
```

hrmcom.txt

```
//
////////////////////////////////////
//
////////

////////////////////////////////////
//
////////
//
//      MONITOR CONNECTION METHODS
//      -----
//
#define          HRMCOM_CONNECTION_UPLINK          0
#define          HRMCOM_CONNECTION_IR              1
//
//      Polar UpLink technology can be used only for transferring
settings from computer
//      to Polar S-series HR monitor (one-way). "Read" functions can b
e
called with
//      HRMCOM_CONNECTION_UPLINK as connection method, but method is
automatically
//      changed to HRMCOM_CONNECTION_IR. Infrared conenction is
automatically two-way,
//      this means all the settings etc. can be read and written.
//
//      When infrared is used for writing or reading data to/from HR
monitor, the communication
//      have to be started by using function
fnHRMCom_StartIRCommunication. After calling this
//      function, all the other reading and writing functions can be
used normally. To end
//      infrared communication, call function
fnHRMCom_EndIRCommunication.
//
////////////////////////////////////
//
////////

////////////////////////////////////
//
////////
//
//      POLAR UPLINK WAVE FILES
//
//      By default wave file (random name to Temp folder) will be
created, played and deleted.
//      Wave file will be created automatically to Temp folder defined
```

hrmcom.txt

at the Windows system.

// The new wave file will be automatically named as HRMxxx.WAV,
where xxx is a random number.

// The playing of wav file do not allow cancelling.

//
////////////////////////////////////
//
/////////
□

////////////////////////////////////
////////////////////////////////////
//
// LIBRARY VERSION DATA
//
////////////////////////////////////
////////////////////////////////////

// Get hrmcom library file version
__declspec (dllexport) int CALLBACK fnHRMCom_GetLibraryVersion
(void);

// Version 1.00 will be returned as 100

////////////////////////////////////
////////////////////////////////////
//
// GENERAL SETTINGS DATA
//
////////////////////////////////////
////////////////////////////////////

// The following data structure will be used with the most of the
functions to give general information
// about communication, for example are we using Polar UpLink or
Infrared connection.

```
typedef struct
{
    int iSize; // Structure
    size for version control // Get

    using sizeof (STRUCTURE)

    int iConnection; // Connection method:
```

```

                                hrmcom.txt
HRMCOM_CONNECTION_UPLINK or HRMCOM_CONNECTION_IR                                // NOT
E:
Polar UpLink connection can be used only for writing information to HR
monitor:

    int            iMonitorID;                                // Unique
monitor ID, 0 = message to all monitors
                                                                //
Monitor will accept the messages if monitor id to send is same as
already
                                                                // set
by User settings or if message was meant for all monitors available.
                                                                // Oth
er
ID numbers used mainly with IR communication

    TCHAR    szWaveFile[MAX_PATH]; // Wave file name, use NULL to
create random file name
                                                                // to
Temp folder (MAX_PATH = 260)
                                                                //
EXCLUSION FLAGS
                                                                //
-----
    BOOL    bNoCreateWave; // Don't create wave file at
all, this allows testing of values in the data structure
    BOOL    bNoPlayWave; // Don't play created Polar
UpLink WAV file
    BOOL    bNoDeleteWave; // Don't delete created Polar
UpLink WAV file after it have been played
                                                                // DAT
A
FILE MANAGEMENT
                                                                //
-----
    BOOL    bLoadFromDataFile; // Load information from binar
y
data file, file name have to be at szWaveFile
                                                                // If
trying to load the data file with not the same data as data structure
                                                                //
specified in call, all the calling functions will return FALSE
                                                                // Whe

```

hrmcom.txt

n
file will be loaded, other actions (create, play, delete wave) are not
done.

loaded file includes incorrect data, default values will be set // If
automatically.

BOOL bSaveAsDataFile; // Save information to binary
data file, file name have to be at szWaveFile

// Whe
n
file will be saved, other actions (create, play, delete wave) are not
done.

CONNECTION DIALOG //

----- //

BOOL bConnectionDlg; // Usage of connection dialog

to
user

HWND hOwnerWnd; // Owner window handle

to connection dialog

connection dialog has been selected to be shown, owner window // If

handle have to be specified. If not, dialog won't be shown and
connection fails. //

connection dialog is not in use, this parameter is ingnored. // If

TCHAR szDlgMsg[50]; // Connection dialog message t
o
user, max 50 characters

message text is not specified, default English texts will be used // If

connection dialog is not in use, this parameter is ingnored. // If

// MIS

C
PARAMETERS

----- //

BOOL bFixErrors; // Errors in settings
can be fixed automatically and error messages

hrmcom.txt

// are

not returned in normal cases.

```
        int          iParam;                // Parameter
reserved for future usage, use zero
        long         lParam;                // Parameter reserved
for future usage, use zero

} POLAR_SSET_GENERAL;
```

```
////////////////////////////////////
////////////////////////////////////
//
//      WATCH SETTINGS
//
////////////////////////////////////
////////////////////////////////////
```

```
//      All Polar S-series HR monitors do have two independent time
zones. The active time zone
//      can be selected with iActiveTime.
```

```
typedef struct
{
        int          iSize;                // Structure
size for version control
                                           // Get
```

using sizeof (STRUCTURE)

```
        int          iTime1;                // Time in
seconds from midnight (0:00:00), max 23:59:59 = 86399 sec
                                           // If
iTime1 = -1, current system time is automatically set to iTime1
        int          iTime2;                // Time in
seconds from midnight (0:00:00), max 23:59:59 = 86399 sec
                                           // If
iTime2 = -1, current system time is automatically set to iTime2
                                           // Onl
```

y
full hours and minutes are valid, seconds will be set to zero

```
        int          iTime1HourMode;        // 0 = 24h mode, 1 = 1
2h
mode
```

```

                                hrmcom.txt
int          iTime2HourMode;          // 0 = 24h mode, 1 = 1
2h
mode
int          iActiveTime;              // 0 = time1 active, 1
=
time2 active
int          iDate;                    // Date in
format yyymmdd, Jan 1 2000 - Dec 31 2099
// If
iDate = -1, current system date is automatically set to iDate
BOOL        bAlarmEnabled;            // FALSE = off, TRUE = on
int          iAlarmTime;               // Time in
seconds from midnight (0:00:00), max 23:59:59 = 86399 sec
// Onl
y
full hours and minutes are valid, seconds will be set to zero

} POLAR_SSET_WATCH;

__declspec (dllexport) void CALLBACK fnHRMCom_ResetWatchSettings
(POLAR_SSET_WATCH*);
__declspec (dllexport) BOOL CALLBACK fnHRMCom_SendWatchSettings
(POLAR_SSET_GENERAL*, POLAR_SSET_WATCH*);
__declspec (dllexport) BOOL CALLBACK fnHRMCom_ReadWatchSettings
(POLAR_SSET_GENERAL*, POLAR_SSET_WATCH*);

// NOTE: Education HR monitors E200 and E600 do have only one time
(iTime1) and no alarm available.
// Set iTime2 to the same as iTime1, hour mode should be the same for
both times.
// Alarm time should be zero and alarm should be not enabled.

□

////////////////////////////////////
////////////////////////////////////
//
// EXERCISE SET
//
////////////////////////////////////
////////////////////////////////////

// Exercise Set information will be send to monitor one set at a time.
// Exercise Set can be set as an active set to monitor (i.e. set will
be
shown
// as the first set when next time starting exercise).

```

hrmcom.txt

```
typedef struct
{
    int          iSize;                // Structure
    size for version control           // Get

using sizeof (STRUCTURE)

    int          iSetNumber;          // Exercise se
t
number: 0, 1, 2, 3, 4, 5.             // Set

number 0 can be used only for setting "Basic Set" active.           // Set
s
2 - 5 are not available for all monitors (see Polar S-series Monitor
Features).

    BOOL        bActiveSet;           // Will this set to be
set as an active set in monitor? TRUE/FALSE
    TCHAR       szName[8];            // Exercise set name
(see Polar Character Set)           //
String can be checked using function fnHRMCom_CheckPolarCharString // Max
number of characters is 7 + ending zero //
"BasicUse" name is permanent for S610, S710 and S810, it can't be
modified.

    // Timers
    BOOL        bTimer1Enabled;        // Timer 1 enabled TRUE/FALSE
    int         iTimer1;               // Timer 1 val
ue
in seconds, max 99 min 59 sec (= 5999 sec)
    BOOL        bTimer2Enabled;        // Timer 2 enabled TRUE/FALSE
    int         iTimer2;               // Timer 2 val
ue
in seconds, max 99 min 59 sec
er                                     // Tim
2 used as interval timer, if intervals enabled.
    BOOL        bTimer3Enabled;        // Timer 3 enabled TRUE/FALSE
    int         iTimer3;               // Timer 3 val
ue
```

hrmcom.txt

in seconds, max 99 min 59 sec

```
// HR Limits
BOOL    bHRLimit1Enabled;    // HR Limits 1 enabled
int      iHRLimit1Upper;    // HR Limit 1 upper
value 30 - 240 bpm
int      iHRLimit1Lower;    // HR Limit 1 lower
value 30 - 240 bpm (must be less than upper limit)
BOOL    bHRLimit2Enabled;    // HR Limits 2 enabled
int      iHRLimit2Upper;    // HR Limit 2 upper
value 30 - 240 bpm
int      iHRLimit2Lower;    // HR Limit 2 lower
value 30 - 240 bpm (must be less than upper limit)
BOOL    bHRLimit3Enabled;    // HR Limits 3 enabled
int      iHRLimit3Upper;    // HR Limit 3 upper
value 30 - 240 bpm
int      iHRLimit3Lower;    // HR Limit 3 lower
value 30 - 240 bpm (must be less than upper limit)

BOOL    bMaxHRInUse;    // Are HR limit values in
percentage of maximum HR given in iMaxHR variable?

// If
TRUE, all HR limit values are used as percentage values (50 - 100%)
int      iMaxHR;    // Maximum HR
value to be used for calculation of HR limit values.

// HR
value in bpm, 100 - 240 bpm

// Intervals
BOOL    bIntervalsEnabled;    // TRUE/FALSE
int      iIntervalType;    // 0 = manual, 1 = timer
er
(use Timer2), 2 = HR,

// 3 =
distance (distance only with cycling models)
int      iIntervalCount;    // The number of
intervals, 0 - 30 (0 = unlimited)
int      iIntervalEndHR;    // Interval ending HR
bpm 10 - 240 bpm
int      iIntervalDistKm;    // The distance of
interval in 0.1 km (max 99.9 km)

// If
monitor does not support cycling features, this value is ignored
int      iIntervalDistMiles;    // The distance of
interval in 0.1 miles (max 99.9 miles)

// If
monitor does not support cycling features, this value is ignored
```

hrmcom.txt

```

// If
both distance values are specified, km value takes precedence.

// Recovery
BOOL      bRecoveryEnabled;      // TRUE/FALSE
int       iRecoveryType;         // 0 = timer recovery,
1
= HR recovery
// 2 =
distance recovery (distance only with cycling models)
int       iRecoveryTime;         // Recovery time in
seconds, max 99 min 59 sec (max 5999 sec)
int       iRecoveryHR;           // recovery HR value 1
0
- 240 bpm
int       iRecoveryDistKm;       // The distance of
recovery in 0.1 km (max 99.9 km)
// If
monitor does not support cycling features, this value is ignored
int       iRecoveryDistMiles;    // The distance of
recovery in 0.1 miles (max 99.9 miles)
// If
monitor does not support cycling features, this value is ignored
// If
both distance values are specified, km value takes precedence.

} POLAR_SSET_EXERCISESET;

__declspec (dllexport) void CALLBACK fnHRMCom_ResetExerciseSet
(int, POLAR_SSET_EXERCISESET*, int, int); // ..., iExerciseType,

iMonitor
__declspec (dllexport) BOOL CALLBACK fnHRMCom_SendExerciseSet
(POLAR_SSET_GENERAL*, POLAR_SSET_EXERCISESET*);
__declspec (dllexport) BOOL CALLBACK fnHRMCom_ReadExerciseSet
(POLAR_SSET_GENERAL*, int, POLAR_SSET_EXERCISESET*);

// Integer value at resetting and reading functions include exercise
number (1 - 5)
// If sending was successful, function returns TRUE, otherwise
FALSE

// NOTE: Education HR monitors E200 and E600 do have only one Exercise
Set, iSetNumber should be 1 (one).
// E200 and E600 HR monitors do have the following features:
```

hrmcom.txt

```
//      - Timers 1 and 2
//      - HR Limits 1
//      - Recovery calculation (type timer always)
// Other exercise settings should be set to default values
□

////////////////////////////////////
////////////////////////////////////
//
// USER SETTINGS
//
////////////////////////////////////
////////////////////////////////////

//      User settings include both information about the person and the
//      usage of monitor features.
//      All the settings are not available in all Polar S-series monitors,
//      see more details
//      from monitor specifications. If any data sent to monitor is not
//      supported, it will
//      be ignored automatically.

typedef struct
{
    int iSize; // Structure
    size for version control // Get

using sizeof (STRUCTURE)
    // Information about user
    int iDateOfBirth; // Date of birth in
format yyymmdd, Jan 1 1921 - Dec 31 2020
    int iActivityLevel; // Activity level: 0 =
low, 1 = moderate, 2 = high, 3 = top
    int iMaxHR; // Maximum heart
rate value 100 - 240 bpm
    int iVO2max; // VO2max value
10 - 95 mmol/l/kg
    int iUserSex; // Sex of user
:
0 = male, 1 = female
```

```

                                hrmcom.txt
    int          iWeightKg;          // Weight in
kilograms: 0, 20 - 199 kg
    int          iWeightLbs;         // Weight in
pounds: 0, 44 - 499 lbs
// If
both weight values are specified, kg value takes precedence.
    int          iHeightCm;          // Height in
centimeters, 0, 90 - 211 cm
    int          iHeightFt;          // Height in
feet: 0, 3 - 7 ft
    int          iHeightInches;      // Height in inches: 0
-
11 inches
// If
both height values are specified, cm value takes precedence.
    TCHAR        szName[8];          // User name (see Polar
r
Character Set)
//
String can be checked using function fnHRMCom_CheckPolarCharString
// Max
number of characters is 7 + ending zero
// If
monitor does not support user name, this value is ignored
// Monitor Features and Functions
    int          iMonitorID;          // Monitor ID
number (for example player number) 0 - 99
    BOOL         bOwnCal;              // OwnCal calculation
enabled TRUE/FALSE
    BOOL         bHRMaxP;              // HRmax-p calculation
enabled TRUE/FALSE
    BOOL         bOwnIndex;            // OwnIndex calculation
n
enabled TRUE/FALSE
    BOOL         bAltimeter;           // Altimeter enabled
TRUE/FALSE, available only for S710
    BOOL         bButtonSound;         // Button sounds enabled
TRUE/FALSE
    BOOL         bOptionsLock;         // Options mode lock enabled
TRUE/FALSE
    BOOL         bHelp;                // Feature help functi
on
enabled TRUE/FALSE
    BOOL         bUS_Units;            // Measurement units:

```

```

                                hrmcom.txt
FALSE = EURO units, TRUE = US units
    int          iSamplingRate;          // 0 = 5s, 1 = 15s, 2
=
60s, 3 = R-R intervals
//
Sampling rate selection is available only with S610, S710 and S810
// R-R
intervals recording is available only with S810
//
Monitor S210 do not have sampling rate selection
//
Monitors S410 and S510 have always dynamic sampling rate
    int          iHeartTouch;          // Usage of Wireless
Button trigger (heart touch feature)
// 0 =
normal, 1 = lap, 2 = change display and limits
//
Wireless button action selection is available with S610, S710 and S810
    int          iRLXBaseLine;          // Relaxation base line
e
only for S810, 4 - 150 mseconds
    BOOL        bOnlineRecording;          // Online recording enabled
TRUE/FALSE, S810 only

} POLAR_SSET_USER;

__declspec (dllexport) void CALLBACK fnHRMCom_ResetUserSettings
(POLAR_SSET_USER *);
__declspec (dllexport) BOOL CALLBACK fnHRMCom_SendUserSettings
(POLAR_SSET_GENERAL*, POLAR_SSET_USER*);
__declspec (dllexport) BOOL CALLBACK fnHRMCom_ReadUserSettings
(POLAR_SSET_GENERAL*, POLAR_SSET_USER*);

//      If sending was succesfull, function returns TRUE, otherwise
FALSE

// NOTE: Education HR monitors E200 and E600 do have only the followin
g
features:
// Options Lock, User Name, Monitor ID, Sampling Rate (E600 only)

////////////////////////////////////
//
////////
//
// VO2max and HRmax-p values are used in OwnCal calories calculation

```

hrmcom.txt

// and those values can be updated as follows:

//
// PC using UpLink/IR ----->
// OwnIndex from FitTest --> UserSet in Monitor ---> OwnCal

calculation in monitor

// Manually set ----->

//
////////////////////////////////////
//
////////

□

////////////////////////////////////
////////////////////////////////////
//
// REMINDER
//
////////////////////////////////////
////////////////////////////////////

// Reminders are available with S610, S710 and S810 heart rate
monitors. There are
// 7 reminder "slots" available in each HR monitor and those can
be
modified only by using
// computer. Each reminder can be individually set to be activate
d
at selected date & time.
// One reminder at time can be sent to HR monitor, select reminde
r
"slot" to be updated by iNumber.
// Reminder can be repeated automatically hourly, daily, monthly,
weekly, monthly and yearly.
// An exercise (ExeSet / ExeProfile) can be set to be active afte
r
reminder has alarmed.

typedef struct
{
 int iSize; // Structure
size for version control // Get

using sizeof (STRUCTURE)

 int iNumber; // Number of

hrmcom.txt

```
reminder, 0 - 6
    BOOL    bActive;                // Reminder activated
TRUE/FALSE
    int      iDate;                // Date of
reminder in format yyyymmdd, Jan 1 2000 - Dec 31 2020
    int      iTime;                // Time in
seconds from midnight (0:00:00), max 23:59:59 = 86399 sec
                                                // Onl
y
full hours and minutes are valid, seconds will be set to zero
    int      iRepeat;              // Repetition
of
reminder. 0 = Off, 1 = Hourly,
                                                // 2 =
Daily, 3 = Weekly, 4 = Monthly, 5 = Yearly
    int      iExercise;            // S810:
Exercise Profile to be set as default profile after reminder alarm
                                                // 0 =
Off, 1 = BasicUse, 2 - 8 Profile Number (remember to update also
exercise profiles)
                                                // S61
0
& S710: Exercise Set to be set as default profile after reminder alarm
                                                // 0 =
Off, 1 = BasicUse, 2 - 7 ExeSet Number (remember to update also
exercise sets)
    TCHAR    szText[8];            // Reminder Text (see
Polar Character Set)
                                                //
String can be checked using function fnHRMCom_CheckPolarCharString
                                                // Max
number of characters is 7 + ending zero

} POLAR_SSET_REMINDER;

__declspec (dllexport) void CALLBACK fnHRMCom_ResetReminder
(int, POLAR_SSET_REMINDER*);
__declspec (dllexport) BOOL CALLBACK fnHRMCom_SendReminder
(POLAR_SSET_GENERAL*, POLAR_SSET_REMINDER*);
__declspec (dllexport) BOOL CALLBACK fnHRMCom_ReadReminder
(POLAR_SSET_GENERAL*, int, POLAR_SSET_REMINDER*);
```

hrmcom.txt

```
//      Integer value at resetting and reading functions include
remainder number (0 - 6)
//      If sending was succesfull, function returns TRUE, otherwise
FALSE

□

////////////////////////////////////
////////////////////////////////////
//
// BIKES
//
////////////////////////////////////
////////////////////////////////////

//      Bike information is available only with S510 and S710. Bike's
power settings will
//      be transferred only to Polar S710 HR monitor.

typedef struct
{
    // Bike Information
    TCHAR    szBikeID[5];          // Bike ID (name) (see Polar
Character Set)

//
String can be checked using function fnHRMCom_CheckPolarCharString
// Max
number of characters is 4 + ending zero
int          iWheelSize;          // Wheel size
in
millimeters (1000 - 3000 mm)
    BOOL      bAutoStart;          // Is autostart featur
e
in use TRUE/FALSE
    BOOL      bSensorSpeed;          // Speed sensor in use
TRUE/FALSE. This flag is not in use, speed sensor is always in use.
    BOOL      bSensorCadence;        // Cadence sensor in use
TRUE/FALSE
    BOOL      bSensorPower;          // Power sensor in use
TRUE/FALSE
// Pow
er
sensor is available only with Polar S710 HR monitor
// If
monitor does not support power sensor, this value is ignored
```

```

                                hrmcom.txt
// Power Sensor Settings
int          iChainMass;           // Weight of
chain in grams (200 - 400 g)
int          iChainLength;        // Length of chain in
mm
(1000 - 2000 mm)
int          iChainWank;          // The length
of
vibrating part (span) chain in mm (300 - 600 mm)

} POLAR_BIKE_INFO;

typedef struct
{
    int          iSize;           // Structure
size for version control
                                // Get

using sizeof (STRUCTURE)

    int          iBikeInUse;      // Which bike
has been selected to be in use right now?
                                // 0 =

Bike1, 1 = Bike2, 2 = None (no cycling features in use)
                                // If
Bike1 is in use, some cycling sensor (speed, cadence or
                                //
power) have to be in use.
    POLAR_BIKE_INFO Bike[2];

} POLAR_SSET_BIKES;

__declspec (dllexport) void CALLBACK fnHRMCom_ResetBikeSettings
(POLAR_SSET_BIKES*);
__declspec (dllexport) BOOL CALLBACK fnHRMCom_SendBikesSettings
(POLAR_SSET_GENERAL*, POLAR_SSET_BIKES*);
__declspec (dllexport) BOOL CALLBACK fnHRMCom_ReadBikesSettings
(POLAR_SSET_GENERAL*, POLAR_SSET_BIKES*);

//      If sending was succesfull, function returns TRUE, otherwise
FALSE

□

////////////////////////////////////
////////////////////////////////////
//

```

hrmcom.txt

```
// EXERCISE PROFILE
```

```
//
```

```
////////////////////////////////////
```

```
////////////////////////////////////
```

```
typedef struct
```

```
{
```

```
    BOOL    bPhaseEnabled;          // Has exercise phase been
enabled? TRUE/FALSE
```

```
    // HR Limits
```

```
    BOOL    bHRLimitEnabled;        // HR Limits enabled
```

```
    int      iHRLimitUpper;          // HR Limit upper valu
```

```
e
```

```
30 - 240 bpm
```

```
    int      iHRLimitLower;          // HR Limit lower valu
```

```
e
```

```
30 - 240 bpm (must be less than upper limit)
```

```
    // Interval period
```

```
    BOOL    bIntervalsEnabled;      // Is entire work period
enabled? TRUE/FALSE
```

```
    int      iIntervalType;          // 0 = manual, 1 =
timer, 2 = End HR
```

```
    int      iIntervalCount;         // The number of
intervals, 0 - 30 (0 = unlimited)
```

```
    int      iIntervalTimer;         // Timer value in
seconds, max 99 min 59 sec (= 5999 sec)
```

```
    int      iIntervalEndHR;         // Interval ending HR
bpm 10 - 240 bpm
```

```
    // Recovery period
```

```
    BOOL    bRecoveryEnabled;        // Is entire recovery period
enabled ? TRUE/FALSE
```

```
    int      iRecoveryType;          // 0 = timer recovery,
1
```

```
= HR recovery
```

```
    int      iRecoveryTime;          // Recovery time in
seconds, max 99 min 59 sec (max 5999 sec)
```

```
    int      iRecoveryHR;            // recovery HR value 1
0
```

```
- 240 bpm
```

```
} POLAR_EXEPHASE;
```

```
typedef struct
```

```
{
```

```
    int      iSize;                  // Structure
```

hrmcom.txt

size for version control

// Get

using sizeof (STRUCTURE)

```
        int            iNumber;                // Exercise
profile number, 1 - 7
        BOOL    bActiveProfile;                // Will this profile to be set
as an active set in monitor? TRUE/FALSE
        TCHAR    szName[8];                    // Exercise profile na
me
(see Polar Character Set)
//
String can be checked using function fnHRMCom_CheckPolarCharString
// Max
number of characters is 7 + ending zero
```

```
        BOOL    bMaxHRInUse;                    // Are HR limit values in
percentage of maximum HR given in iMaxHR variable?
// If
TRUE, all HR limit values are used as percentage values (50 - 100%)
        int            iMaxHR;                // Maximum HR
value to be used for calculation of HR limit values.
// HR
value in bpm, 100 - 240 bpm
```

```
        POLAR_EXEPHASE    Phase[6];            // One exercise profile includ
es
6 exercise phases
// Eac
h
phase should be defined as POLAR_EXEPHASE structure
```

```
} POLAR_SSET_EXERCISEPROFILE;
```

```
__declspec (dllexport) BOOL CALLBACK fnHRMCom_ResetExerciseProfile
(int, POLAR_SSET_EXERCISEPROFILE*, int);
```

```
__declspec (dllexport) BOOL CALLBACK fnHRMCom_SendExerciseProfile
(POLAR_SSET_GENERAL*, POLAR_SSET_EXERCISEPROFILE*);
```

```
__declspec (dllexport) BOOL CALLBACK fnHRMCom_ReadExerciseProfile
(POLAR_SSET_GENERAL*, int, POLAR_SSET_EXERCISEPROFILE*);
```

```
//      NOTE: Exercise profiles are available only with Polar S810 HR
monitor.
```

```
//      Integer value at resetting and reading functions include Exe
Profile number (1 - 7)
```

```
//      If sending was succesfull, function returns TRUE, otherwise
FALSE
```

```
□
```

```
//////////
//////////
//
// MONITOR BITMAP LOGO
//
//////////
//////////
```

```
//      NOTE: Monitor bitmap logos are available with Polar S610, S710
,
S810, E200 and E600 HR monitors.
```

```
__declspec (dllexport) BOOL CALLBACK fnHRMCom_SendBitmap
```

```
(POLAR_SSET_GENERAL*, int*);
```

```
__declspec (dllexport) BOOL CALLBACK fnHRMCom_ReadBitmap
```

```
(POLAR_SSET_GENERAL*, int*);
```

```
//      If sending was succesfull, function returns TRUE, otherwise
FALSE
```

```
// Example:      Each pixel column in one integer value => int
iBitmapPixelCol[47];
//                  First pixel in the bottom is 2^0, second 2^1,
third 2^2, etc.
//                  If three pixels in bottom are ON =>
iBitmapPixelCol[iColumn] = 7 (1+2+4)
//                  If entire column is ON =>
iBitmapPixelCol[iColumn] = 255 (1+2+4+8+16+32+64+128)
//                  Send to monitor fnHRMCom_SendBitmap
(&iBitmapPixelCol[0]);
```

```
□
```

```
//////////
//////////
//
// SET HR MONITOR TO WATCH MODE
//
//////////
//////////
```

hrmcom.txt

```
//      Sets monitor to watch mode, monitor do not accept other
//      messages, until it has been switched back to Connect mode.
```

```
__declspec (dllexport) BOOL CALLBACK fnHRMCom_SendMonitorToWatchMode
(POLAR_SSET_GENERAL*);
```

```
//      If sending was succesfull, function returns TRUE, otherwise
FALSE
```

□

```
////////////////////////////////////
////////////////////////////////////
//
//  FACTORY DEFAULTS
//
////////////////////////////////////
////////////////////////////////////
```

```
//      Sets monitor factory defaults, resets all monitor data includi
ng
EEPROM memory. Use very carefully!!!
//      Setting factory defaults is not meant for normal software usag
e,
only for service software products.
//      When settings factory defaults, confirmation of the operation
should be asked always.
```

```
__declspec (dllexport) BOOL CALLBACK fnHRMCom_SendFactoryDefaultComman
d
(POLAR_SSET_GENERAL*);
```

□

```
////////////////////////////////////
////////////////////////////////////
//
//  DELETING EXERCISE FILES FROM HR MONITOR
//
////////////////////////////////////
////////////////////////////////////
```

```
//      Exercise files can be deleted by using the following function
call. Files can be deleted
//      from Polar S610, S710, S810 and E600 HR monitors.
//
//      NOTE: The first version of Polar S610 (DataVersion=1) can't
handle deleting one exercise
```

```

                                hrmcom.txt
//      file correctly, all exercise files can still be deleted.
//
//      NOTE: The first version of Polar S810 (DataVersion=3) can't
handle deleting all exercise
//      files correctly, one exercise file can still be deleted.
//
//      Check monitor type and data version before sending file delete
message to monitor!
//
//      Give exercise number as int parameter (0 - n), if all files
should be deleted, use int parameter -1

__declspec (dllexport) BOOL CALLBACK fnHRMCom_SendFileDeleteCommand
(POLAR_SSET_GENERAL*, int);

□

////////////////////////////////////
////////////////////////////////////
//
// INFRARED COMMUNICATION FUNCTIONS
//
////////////////////////////////////
////////////////////////////////////

////////////////////////////////////
//
// MONITOR INFO
//
////////////////////////////////////

typedef struct
{
    int                iSize;                                //
    Structure size for version control

    int                iMonitorInUse;                        // HR monitor
in
use: HRM_S610, HRM_S710, HRM_S810 or HRM_E600
    int                iDataVersion;                        // HR monitor
data version

    int                iTotFiles;                            // Total count
of all files inside HR monitor
    int                iFreeMemoryInBytes;                  // Free memory

```

```

                                hrmcom.txt
inside HR monitor (in bytes)
    int                iTotalMemoryInBytes;    // Total memory inside
HR monitor (in bytes)

    BOOL    bLowBattery;                        // Low battery indicat
or
TRUE / FALSE

} POLAR_SSET_MONITORINFO;

__declspec (dllexport) BOOL CALLBACK fnHRMCom_ReadMonitorInfo
(POLAR_SSET_GENERAL*, POLAR_SSET_MONITORINFO*);

□

////////////////////////////////////
//
////////
//
//      fnHRMCom_ResetIRCommunication
//      -----
//      Call in the startup of software to reset all the communication
parameters.
//
//      Parameters:
//      int iParam                Reserved in future usage, use 0 (zero)
.
//
//      Return value:
//      TRUE                    - Resetting made succesfully
//      FALSE                   - Resetting was not made because of
communication is already running.
//
////////////////////////////////////
//
////////
//

__declspec (dllexport) BOOL CALLBACK fnHRMCom_ResetIRCommunication
(int);

□

////////////////////////////////////
//
////////
//

```

```

                                hrmcom.txt
//      fnHRMCom_StartIRCommunication
//      -----
//      Call to open communication port and start IR communication.
//
//      Parameters:
//      int iParam                Parameter for connection settings (when
used multiple params, use OR)
//
//                                See "Infrared communication
parameters"
//
//      LPTSTR tcPort            Communication port name, for example "COM1:" o
r
"COM2:"
//
//                                Remember to use to use colon :
at the end of port name
//
//      Return value:
//      BOOL bStartOK
//      TRUE                    - Starting of communication made
successfully
//      FALSE                    - Problems encountered, check the
following possible errors:
//
//                                * Communication has
already been started and it is running
//                                * Communication port
already reserved for some other device
//                                * Maybe call was made
from 16-bit program. A 32-bit DLL cannot
//                                create an additional
thread when that DLL is being called by
//                                a 16-bit program.
//
////////////////////////////////////
//
/////////

__declspec (dllexport) BOOL CALLBACK fnHRMCom_StartIRCommunication
(int, LPTSTR);

□

// Infrared communication parameters
#define HRMCOM_PARAM_INTERNALIR 1 // Use interna
l
IR port (Win95 only)

```

```

                                hrmcom.txt
#define HRMCOM_PARAM_KEEPCONNECT          2          // Keep monitor
in Connect mode during connection
#define HRMCOM_PARAM_FILTERHRDATA         4          // Filter
averaged HR data (sampling rate 5 or 15 seconds)

// (not yet implemented)
#define HRMCOM_PARAM_DIRECT_USB           16         // Direct USB port usage
#define HRMCOM_PARAM_VIRTUALCOMPORT       32         // Virtual COM port usage
#define HRMCOM_PARAM_DUMPFRAMES           64         // Dump frames to
c:\frames.txt or c:\all.txt text files

// Dumping can be used for data error detection.
#define HRMCOM_PARAM_ONLINE                128        // Online
recording mode (Polar S810 only)

□

////////////////////////////////////
//
////////
//
//      fnHRMCom_EndIRCommunication
//      -----
//      Call to close communication port and end IR communication.
//
//      Parameters:
//      int iParam          - Reserved in future usage, use 0
//      (zero).
//
//      Return value:
//      BOOL bEndOK
//      TRUE                - Ending of communication made
//      successfully
//      FALSE               - Problems with ending of communication
//
////////////////////////////////////
//
////////

__declspec (dllexport) BOOL CALLBACK fnHRMCom_EndIRCommunication
(int);

```

hrmcom.txt

□

```

////////////////////////////////////
//
////////
//
//      Communication Texts
//      -----
//
//      Communication texts are shown with infrared communication
process. By default
//      English texts for buttons and message texts are defined. If
texts need localization,
//      the following functions can be used to change communication
texts before calling
//      communication functions. The text at the end of the following
defines shows the
//      default text string for each text item.
//
#define          HRMCOM_TEXT_CANCEL                      0

// Cancel
#define          HRMCOM_TEXT_RETRY                        1

// Retry
#define          HRMCOM_TEXT_READING                     2

// Reading...
#define          HRMCOM_TEXT_NOANSWER                     3
// No answer from HR Monitor
#define          HRMCOM_TEXT_ERRORS                       4

// Errors with Connection
#define          HRMCOM_TEXT_STARTING                     5
// Starting Connection...
#define          HRMCOM_TEXT_TITLE                        6

// Infrared Connection
#define          HRMCOM_TEXT_WRITING                      7

// Writing...
//
//      To set each communication text, call function
fnHRMCom_SetComText.
//      For example this call will change internal text for informing
user about
//      not getting any answers from HR monitor within answer time:
//

```

```

                                hrmcom.txt
//      fnHRMCom_SetComText (HRMCOM_TEXT_NOANSWER, "Ei vastausta
sykemittarilta");
//
__declspec (dllexport) BOOL CALLBACK fnHRMCom_SetComText
(int, LPTSTR);
//
//      When HRMCOM.DLL is initialized by starting software calling it
,
all the communication
//      texts are resetted automatically. To reset all communication
texts back to English
//      default texts, the following reset function can be used:
//
__declspec (dllexport) void CALLBACK fnHRMCom_ResetComTexts (void)
;
//
//      NOTE: Title text for Polar UpLink Communication is always "Pol
ar
UpLink".
//
////////////////////////////////////
//
/////////

□

////////////////////////////////////
//
/////////
//
//      Reading Exercises Data
//      -----
//
//      Reading exercises data from HR monitor using IR connection can
be done by
//      sending and answering to each communication message and also b
y
calling function,
//      which generates communication dialog and handles all message b
y
itself.
//      Communication port has to be opened before reading exercises
from HR monitor.
//      To read all exercises to memory of DLL, use the following
function:
//
__declspec (dllexport) BOOL CALLBACK fnHRMCom_ReadExercisesData

```

hrmcom.txt

```
(HWND, BOOL);
//
//      Parameters:
//      HWND      hOwnerWnd      - Handle to owner window
//      BOOL      bOneWay        - Flag for one way connection
//      (under construction)
//
//      After all exercises have been read from HR monitor, the basic
//      information
//      about each exercise can be read by using the following function
fnHRMCom_GetExeFileInfo.
//      Structure POLAR_EXERCISEFILE includes basic information about
//      the exercise data
//      file requested.
//
typedef struct
{
    int          iSize;                      //
    Structure size for version control

    int          iTime;                      // Start
    time of exercise in seconds

    int          iDate;                      // Start
    date of exercise in yyymmdd

    int          iDuration;                  //
    Duration of exercise in seconds

    BOOL         bUSTimeMode;                // Usage of 12h time
    mode in exercise

    int          iSamplingRate;              // Sampling rate
    of exercise

    BOOL         bDeleted;                   // Exercise has
    been marked to be deleted

    BOOL         bSpeed;                     // Speed sensor
    data available

    BOOL         bCadence;                   // Cadence
    sensor data available

    BOOL         bAltitude;                  // Altitude
    sensor data available

    BOOL         bPower;                     // Power sensor
    data available
```

hrmcom.txt

```

        BOOL      bInterval;                                // Interval da
ta
available

        TCHAR      szName[9];                                // Exercise
set/profile name used in exercise

// Max number of characters is 8 + ending zero

} POLAR_EXERCISEFILE;
//
__declspec (dllexport) BOOL CALLBACK fnHRMCom_GetExeFileInfo
(int, POLAR_EXERCISEFILE*);
//
//      Parameters:
//      int iExercise                                Parameter for specifyi
ng
exercise of which the information will be retrieved
//      POLAR_EXERCISEFILE* pef*                    Address to exercise file
information data structure
//
//      Before reading detailed exercise information from HRMCOM.DLL's
memory, each
//      exercise file have to be analyzed by using the following
function:
//
__declspec (dllexport) BOOL CALLBACK      fnHRMCom_AnalyzeFile
(int, int);
//
//      Parameters:
//      int iExercise      Parameter for specifying exercise to be analys
ed
//      int iAction        Parameter for specifying the actions t
o
be doen for analysed HR file
//
//                                (when used multiple params, us
e
OR)
//
//                                HRMCOM_PARAM_FILTERHRDATA
= Filter averaged HR data (sampling rate 5 or 15 seconds)
//

Not available yet!
//
//      After the succesfull analyzing, all the exercise information c

```

hrmcom.txt

```
an
be read
//      by using the functions and defines shown in the following
chapters.
//
////////////////////////////////////
//
////////

□

////////////////////////////////////
//
// HRM DATA OUTPUT FUNCTIONS
//
////////////////////////////////////

__declspec (dllexport) int CALLBACK    fnHRMCom_GetRecParam
(int);          // returns recording parameters
__declspec (dllexport) BOOL CALLBACK   fnHRMCom_GetRecFlags
(int);          // returns recording flags

__declspec (dllexport) int CALLBACK    fnHRMCom_GetNbrOfHRMSamples
(void);         // returns nbr. of samples
__declspec (dllexport) int CALLBACK    fnHRMCom_GetHRMSamples
(int, int);     // returns HR/CC samples

__declspec (dllexport) int CALLBACK    fnHRMCom_GetNbrOfIntTimes
(void);         // returns number of lap times
__declspec (dllexport) int CALLBACK    fnHRMCom_GetIntTimeData
(int, int);     // returns lap time data

__declspec (dllexport) BOOL CALLBACK   fnHRMCom_GetNbrOfSwapTimes
(void);         // returns number of HR limit swaps
__declspec (dllexport) int CALLBACK    fnHRMCom_GetLimitSwapData
(int, int);     // returns limit swap data

////////////////////////////////////
//
// HRM DATA FLAGS
//
//      Get these parameters by using function: fnHRMCom_GetRecFlags
//
////////////////////////////////////

#define          FLAG_CYCLO_DATA          3

// TRUE, cycling data
```

```

                                hrmcom.txt
#define          FLAG_3LIMITS_IN_USE                                6

// three HR limits has been used
#define          FLAG_SPEED_DATA                                    8

// file has speed data
#define          FLAG_ALT_DATA                                     9

// file has altitude data
#define          FLAG_CAD_DATA                                    10

// file has cadence data
#define          FLAG_POWER_DATA                                  11

// file has power data
#define          FLAG_INTERVAL_DATA                               12

// file has interval data
#define          FLAG_LAP_DATA                                    13

// file has lap data
#define          FLAG_LIMSWAP_DATA                                14

// file has limit swap data
#define          FLAG_POWER_BALANCE                               18

// file has LR balance data
#define          FLAG_POWER_INDEX                                 19

// file has pedalling index data

////////////////////////////////////
//
// HRM DATA GENERAL RECORDING INFORMATION
//
//      Get these parameters by using function: fnHRMCom_GetRecParam
//
////////////////////////////////////

#define          REC_AM_PM
1                                // 0 = AM, 1 = PM
#define          REC_MONITOR_TYPE                                5

// HR Monitor Type
#define          REC_EURO_US_UNITS                                8

// 0 = Euro, 1 = US

```

```

                                hrmcom.txt
#define          REC_START_DATE                                     9

// Exercise start date in yyyymmdd format
#define          REC_START_TIME                                    10

// Exercise start time hh:mm:ss.s/10 in 1/10 of seconds
#define          REC_REC_LENGTH                                    11

// Duration on exercise (in ms)
#define          REC_SAMPLING_RATE                                12

// Recording rate
#define          REC_UPPER_LIMIT_1                                13

// 0 - 250 bpm
#define          REC_LOWER_LIMIT_1                                14

// 0 - 250 bpm
#define          REC_UPPER_LIMIT_2                                15

// 0 - 250 bpm
#define          REC_LOWER_LIMIT_2                                16

// 0 - 250 bpm
#define          REC_UPPER_LIMIT_3                                17

// 0 - 250 bpm
#define          REC_LOWER_LIMIT_3                                18

// 0 - 250 bpm
#define          REC_ANAEROB_LIMIT                                19

// 0 - 250 bpm
#define          REC_AEROB_LIMIT                                  20

// 0 - 250 bpm
#define          REC_TIMER_1
21                                // timer 1 in seconds
#define          REC_TIMER_2
22                                // timer 2 in seconds
#define          REC_TIMER_3
23                                // timer 3 in seconds
#define          REC_MAX_HR
25                                // UpperLimit+1 - 250
#define          REC_REST_HR
26                                // 0 - LowerLimit-1
#define          REC_RR_START_DELAY                                27

```

```

// R-R recording start delay
#define          REC_START_SAMPLE          29

// 0 - 250 bpm

#define          REC_STOP_TIME             30

// hh:mm:ss.s/10 in 1/10 of seconds
#define          REC_STOP_SAMPLE           31

// 0 - 250
#define          REC_STOP_SPEED            32

// stop speed
#define          REC_STOP_CAD              33

// stop cadence
#define          REC_STOP_ALT              34

// stop altitude
#define          REC_MIN_HRATE             35

// lowest heart rate
#define          REC_AVE_HRATE             36

// average heart rate
#define          REC_MAX_HRATE             37

// highest heart rate

#define          REC_TRIP_DIST_STOP        38

// trip distance at stop
#define          REC_TRIP_CLIMB_STOP       39

// trip climb at stop
#define          REC_TOT_TIME_STOP         40

// total time at stop
#define          REC_AVG_ALT               // average altitude
41
#define          REC_MAX_ALT               // maximum altitude
42
#define          REC_AVG_SPEED             43

// average speed
#define          REC_MAX_SPEED             44

```

hrmcom.txt

```
// maximum speed
#define REC_ODOM_STOP 45

// odometer stop
#define REC_MIN_SPEED 46

// minimum speed

#define REC_RECOVERY_TIME 47
#define REC_RECOVERY_HR 48
#define REC_MAX_POWER 78

// Maximum power in watts
#define REC_AVE_POWER 79

// Average power in watts
#define REC_CALORIES 80

// Calory consumption
#define REC_NBR_OF_LIMITS_IN_USE 83

// Nbr. of HR limits in use

////////////////////////////////////
//
// HRM DATA SAMPLE TYPES
//
// Before getting measured values (samples), get the number of
samples by using
// function fnHRMCom_GetNbrOfHRMSamples. After this operation,
samples can be get
// by calling function fnHRMCom_GetHRMSamples for example in the
following way:
//
// iTotal = fnHRMCom_GetNbrOfHRMSamples ();
//
// for (i = 0; i < iTotal; i++)
// {
//     iHR[i] = fnHRMCom_GetHRMSamples (CC_HRATE, i);
//     iSpeed[i] = fnHRMCom_GetHRMSamples (CC_SPEED, i);
//     iCad[i] = fnHRMCom_GetHRMSamples (CC_CAD, i);
// }
//
// Speed and altitude values unit depends of recording parameter
REC_EURO_US_UNITS.
// To get the correct units, use for example the following call:
//
```

```

                                hrmcom.txt
//      if (1 == fnHRMCom_GetRecParam (REC_EURO_US_UNITS))
//      {
//          Speed in mph, altitude in feet
//      }
//      else
//      {
//          Speed in km/h, altitude in meters
//      }
//
////////////////////////////////////

#define          CC_HRATE
1                // heart rate values (bpm / msec)
#define          CC_SPEED
2                // speed values (10 * km/h / 10 * mph)
#define          CC_CAD
3                // cadence values (rpm)
#define          CC_ALT
4                // altitude values (m / ft)
#define          CC_POWER
5                // power values (Watts)
#define          CC_POWER_BALANCE                                     6

// power LR Balance (left%)
#define          CC_POWER_INDEX                                     7

// power pedalling index (%)

////////////////////////////////////
//
// LAP TIME DATA INFORMATION
//
//      Before getting lap time data, get the number of laps by using
//      function fnHRMCom_GetNbrOfIntTimes. After this operation, lap
//      information can be get by calling function
fnHRMCom_GetIntTimeData
//      for example in the following way:
//
//      iTotal = fnHRMCom_GetNbrOfIntTimes ();
//
//      for (i = 0; i < iTotal; i++)
//      {
//          iTime   = fnHRMCom_GetIntTimeData (i, INT_INT_TIME);
//          iHR     = fnHRMCom_GetIntTimeData (i, INT_SAMPLE);
//          iSpeed  = fnHRMCom_GetIntTimeData (i, INT_SPEED);
//      }
//

```

hrmcom.txt

```

////////////////////////////////////

```

```

#define          INT_INT_TIME                                601

```

```

// Lap time in 1/10 seconds
#define          INT_LAP_INTRVAL                            603

```

```

// Lap type: 0 = normal lap, 1 = interval
#define          INT_LAP_DISTANCE                            604

```

```

// Lap distance in meters / yards
#define          INT_SAMPLE
607
#define          INT_MIN_SAMPLE                            608
// Momentary HR, 0 - 250 bpm

```

```

// Lap's min HR, 0 - 250 bpm
#define          INT_AVE_SAMPLE                            609

```

```

// Lap's avg HR, 0 - 250 bpm
#define          INT_MAX_SAMPLE                            610

```

```

// Lap's max HR, 0 - 250 bpm
#define          INT_SPEED
611
#define          INT_AVG_SPEED                            612
// Momentary speed, 10 * km/h or mph

```

```

// Average speed, 10 * km/h or mph
#define          INT_CADENCE
613
#define          INT_AVG_CADENCE                            614
// Momentary cadence, 0 - 180 rpm

```

```

// Average cadence, 0 - 180 rpm
#define          INT_ALTITUDE                            615

```

```

// Momentary altitude, (-1000 - 2047) * 10 m / ft
#define          INT_AVG_ALTITUDE                            616

```

```

// Average altitude, (-1000 - 2047) * 10 m / ft
#define          INT_POWER
617
#define          INT_MAX_POWER                            618
// Momentary power, 0 - 2000 Watts

```

```

// Maximum power, 0 - 2000 Watts
#define          INT_AVE_POWER                            619

```

```

// Average power, 0 - 2000 Watts
#define          INT_TEMP
621
// Momentary temperature, 10 * -100 -

```

hrmcom.txt

```
+100 'C or 'F (only from S710)
#define          INT_DIST_REC                      624

// Distance recovery, 10 * km or miles
#define          INT_RECOVERY                      625

// Recovery calculation, 0 = No recovery, 1 = Time Recovery, 2 = +HR
Recovery, 3 = -HR Recovery, 4 = Distance Recovery
#define          INT_HR_REC
626              // HR recovery value, 0 - 3599 seconds
#define          INT_TIME_REC                      627

// Time recovery value, 0 - 240 bpm
#define          INT_LAP_ASCENT                    636

// Lap ascent, trip up,  m / feet

////////////////////////////////////
//
////////
//
// HR LIMIT SWAPS, indexes for 'fnHRMCom_GetLimitSwapData'
//
////////////////////////////////////
//
////////

#define          LIM_SWAP_TIME                      900

// HR limit swap time
#define          LIM_SWAP_CODE                      901

// HR limit swap code

////////////////////////////////////
////////////////////////////////////
//
// ONLINE RECORDING
//
////////////////////////////////////
////////////////////////////////////

//      Online recording is available only with Polar S810 HR monitor.

To start online recording,
//      function fnHRMCom_StartIRCommunication have to be called with
parameter HRMCOM_PARAM_ONLINE.
//      Function 'fnHRMCom_GetOnlineData' returns online data samples
```

hrmcom.txt

```
received from the S810
//      HR monitor or ONLINE_BUFF_EMPTY if there aren't any new online
samples in the buffer.
//
//      iData = fnHRMCom_GetOnlineData (iParam);
//
//      Parameter 'iParam' is 32-bit integer and it is reserved for
future use and it should
//      be 0 (zero) now. Return value is 32-bit integer and it is R-R
value in milliseconds
//      or ONLINE_BUFF_EMPTY if there aren't any new samples in buffer
.
//
////////////////////////////////////
#define      ONLINE_BUFF_EMPTY      -1

// online buffer is empty

__declspec (dllexport) int  CALLBACK      fnHRMCom_GetOnlineData
(int);      // return online data samples
```